

# chattanooga®

## INTELECT® RPW LITE

### EASY ACCESS TO RADIAL PRESSURE WAVE THERAPY

#### EFFICIENCY

Treatment outcomes are determined by the level of energy released by the high quality transmitters. The Intellect® RPW Lite is delivered with the D20-S transmitter. D20-T, Ro40, C15 and DI15 are available as an option.

#### PORTABILITY

Compact and easy to handle.

#### INDICATIONS

- ACHILLES TENDINOPATHY (1)
- PLANTAR FASCIITIS (3)
- TRIGGER FINGER (4)
- TRIGGER POINTS (5)
- PAIN AND INFLAMMATION MANAGEMENT IN ORTHOPAEDIC CONDITIONS (6)

#### TECHNICAL INFORMATION

30680-230V: 230VAC 50-60Hz 700VA

30680-115V: 115VAC 60Hz 700VA

30680-100V: 100VAC 60Hz 700VA

DIMENSIONS: 289x238x310mm

WEIGHT: 9.8Kg

ENERGY: 6 Energy Levels: 6-18Hz/max 2.7 bar

CONFORMITY: this device complies with the applicable standards EN/IEC 60601-1, CAN/CSA-C22.2 No.601.1, UL Std. No 60601-1.

CERTIFICATION: CE 0197

CLASSIFICATION ACCORDING TO MDD: Class IIa device 93/42 60601-1



1. Gerdesmeyer L. Current evidence of extracorporeal shock wave therapy in chronic Achilles tendinopathy. Int J Surg. 2015 \*

2. Ogden JA. Shockwave therapy for chronic proximal plantar fasciitis Foot Ankle Int. 2002 Apr

3. Yildirim P.Extracorporeal shock wave therapy versus corticosteroid injection in the treatment of trigger finger: J Hand Surg Eur Vol. 2016 Nov

4. Nedelka T.Mechano-transduction effect of shockwaves in the treatment of lumbar facet joint pain:Neuro Endocrinol Lett. 2014.

5. Schmitz C. Efficacy and safety of extracorporeal shock wave therapy for orthopedic conditions.Br Med BulL. 2015.



# INTELECT® RPW LITE

## ORDERING INFORMATION

PART NUMBER	DESCRIPTION
30680-230V	Intellect® RPW Lite 230V set including standard accessories
30680-115V	Intellect RPW Lite 115V unit set including standard accessories
30680-100V	Intellect RPW Lite 100V unit set including standard accessories

## STANDARD ACCESSORIES

PART NUMBER	DESCRIPTION
30885-1	Intellect RPW Lite SPARROW™ handpiece kit grey, including 1x projectile; 1x guiding tube; 1x wrench; 1x allen key
	1x Power cord with EU plug for 30680-230V, 1x Power cord with North American plug for 30680-115V, 1x Power cord with Japanese plug for 30680-100V

30873	Transmitter D20-S Black RPW Lite, 20mm
30891	Operating manual and quick start guide

22601	Transmission gel, 250 ml
28739	Guiding tube cleaning brush
30890	Tube set (includes 3 tubes)

## OPTIONAL ACCESSORIES

PART NUMBER	DESCRIPTION
30885-1	Intellect RPW Lite SPARROW handpiece kit grey, including 1x projectile; 1x guiding tube; 1x wrench; 1x allen key
30885-2	Intellect RPW Lite SPARROW handpiece kit green, including 1x projectile; 1x guiding tube; 1x wrench; 1x allen key
29104	SPARROW handpiece revision kit, including 2x projectiles; 2x guiding tubes short; 2x o-ring set; 1x o-ring guide

30873	Transmitter D20-S Black RPW Lite, 20mm
30876	Transmitter D20-T Black RPW Lite, 20mm
30842	Transmitter Ro40 Black RPW Lite, 15mm
30875	Transmitter C15 Black RPW Lite, 15mm
30874	Transmitter DI15 Black RPW Lite, Deep Impact 15mm
30883	Handpiece holder black

30890	Tube set (includes 3 tubes)
22601	Transmission gel, 250 ml
28739	Guiding tube cleaning brush
0.0032.012	Power cord EU
0.0032.012-US	Power cord USA
13-7458-CN	Power cord China/Australia/New Zealand
13448-CH	Power cord Switzerland



DJO® Global | 1a Guildford Business Park | Guildford | Surrey | GU2 8XG | UNITED KINGDOM

[chattanoogarehab.com](http://chattanoogarehab.com)

SPARROW™ is a registered trade mark of Storz medical AG

\*Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Information related to various health, medical, and fitness conditions and their treatment is not meant to be a substitute for the advice provided by a physician or other medical professional.\*

Copyright © [2021] by DJO, LLC All rights reserved. No part of this publication may be reproduced, distributed, or transmitted in any form or by any means.